ZHAO Q, LIU M, ZHANG X, LIN C, ZHANG Q, SHEN Z. Generation of insect-resistant and glyphosate-tolerant rice by introduction of a T-DNA containing two Bt insecticidal genes and an *EPSPS* gene. *Journal of Zhejiang University Science B*. 2015;16(10):824-831. doi:10.1631/jzus.B1500056.

Generation of insect-resistant and glyphosate-tolerant rice by introduction of a T-DNA containing two Bt insecticidal genes and an *EPSPS* gene*

Qi-chao ZHAO, Ming-hong LIU, Xian-wen ZHANG, Chao-yang LIN, Qing ZHANG, Zhi-cheng SHEN^{†‡}

(State Key Laboratory of Rice Biology, Institute of Insect Sciences, Zhejiang University, Hangzhou 310058, China)

[†]E-mail: zcshen@zju.edu.cn

Received Mar. 9, 2015; Revision accepted June 9, 2015; Crosschecked Sept. 13, 2015

Abstract: Insect resistance and glyphosate tolerance have been two of the most important traits in the genetic improvement of various crops. In this study, two *Bacillus thuringiensis* (Bt) insecticidal genes, *Cry1Ac* and *Cry1Ig*, and a modified glyphosate-tolerant 5-enolpyruvylshikimate-3-phosphate synthase (*EPSPS*) gene (*G10*) were combined into a single transferred DNA (T-DNA) fragment and introduced into rice by *Agrobacterium*-mediated transformation. A transgenic line with single-copy T-DNA insertion named GAI-14 was found to be highly resistant to striped stem borer and rice leaf roller, and tolerant to glyphosate. Analysis of T-DNA border sequence suggested that the transgenes were inserted at the chromosome 3 and appeared to have not interrupted any known or putative genes. A field trial observed no significant difference in the basic agronomic traits between GAI-14 and the recipient rice.

Key words: Bt gene stacking, Bt resistance management, *EPSPS*, Transgenic rice

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CDS 278..3814

/gene="cry1Ac"

/note="derived from Bacillus thuringiensis"

/codon_start=1

/product="Cry1Ac"

/protein_id="ADO64599.1"

/translation="MDNNPNINECIPYNCLSNPEVEVLGGERIETGYTPIDISLSLTQ

FLLSEFVPGAGFVLGLVDIIWGIFGPSQWDAFLVQIEQLINQRIEFARNQAISRLEG
LSNLYQIYAESFREWEADPTNPALREEMRIQFNDMNSALTTAIPLFAVQNYQVPLLSV
YVQAANLHLSVLRDVSFVGQRWGFDAATINSRYNDLTRLIGNYTDHAVRWYNTGLERV
WGPDSRDWIRYNQFRRELTTLVLDIVSLFPNYDSRTYPIRTVSQLTREIYTNPVLENF
DGSFRGSAQGIEGSIRSPHLM DILNSITIYTD AHRGEYYWSGHQIMASPVGFSGPEFT
FPLYGTMGNAAPQQRIVAQLGQG VYRTL SSTLYRRPFNIGINNQQLSVLDGTEFAYGT
SSNLPSAVYRKSGTVDSLDEIPPQNNNVPPRQGFSHRLSHVSMFRSGFSNSSVSIIRA
PMFSWIHRSAEFNNIIASDSITQIPAVKGNFLNGSVISGPGFTGGDLVRLNSSGNNI
QNRGYIEVPIHFPSTSTRYRVRVRYASVTPIHLNVNWGNSSIFSNTVPATATSLDNLQ
SSDFGYFESANAFTSSLGNIVGVRNFSGTAGVIIDRFEFIPVTATLEAEYNLERAQKA
VNALFTSTNQLGLKTNVTDYHIDQVSNLVTYLSDEFCLDEKGELSEKVKHAKRLSDER
NLLQDSNFKDINRQPERGWGGSTGITIQQGDDVFKENYVTLSGTFDECYPTYLYQKID
ESKLKAFTRYQLRGYIEDSQDLEIYSIRYNAKHETVNVPGTGSLWPLSAQSPIGKCGE
PNRCAPHLEWNPDLDCSRDGEKCAHHSHHFSLDIDVGCTDLNEDLGVWVIFKIKTQD
GHARLGNLEFLEEKPLVGEALARVKRAEKKWRDKREKLEWETNIVYKEAKESVDALFV
NSQYDQLQADTNIAMIIHAADKRVHSIREAYLPELSVIPGVNAAIFEELEGRIFTAFSL
YDARNVIKDGGFNNGLSRWNVKGHV DVEEQNNQRSVLVVP EWAEVSQEVRVCPGRGY
ILRV TAYKEGYGEGCVTIHEIENNTDELKFSNCVEEEIYPNNTVTCNDYTVNQEEYGG
AYTSRNRGYNEAPSPADYASVYEEKSYTDGRRENPCFNRGYRDTPLPVGCVTKEL
EYFPETDKVWIEIGETEGTFIVDSVELLLMEE"

Cry1lg

source

1..2160

/organism="Bacillus thuringiensis"

/mol_type="genomic DNA"

/strain="ARP166"

/isolation_source="soil"

/db_xref="taxon:1428"

CDS

1..2160

/codon_start=1

/transl_table=11

/product="pesticidal protein"

/protein_id="AGU13866.1"

/translation="MKLKNQDKHQSFSSNAKVDKISTDSLKNETDIELQNINHEDCLK
MSEYENVEPFVSVSTIQTGIGIAGKILGNLGVPFAGQVASLYSFILGELWPKGKSQWE
IFMEHVEELINQKISTYARNKALADLKGLGDALAVYHESLESWIKNRNNTTRTRSVVKS
QYITLELMFVQSLPSFAVSGEEVPLLPIYAQAANLHLLLLRDASIFGKEWGLSDSEIS
TFYNRQVERTSDYSDHCTKWFDTGLNRLKGSNAEIWVKYNQFRRDMTLMVLDLVALFQ
SYDTHMYPIKTTAQLTREVTNAIGTVHPHPSFASTTWYNNNAPSFAIEAAVIRSPH
LLDFLEQVTIYSLLSRWSNTQYMNMWGGHKLEFRTIGGT LNTSTQGSTNTSINPVTLP
FTSRDIYRTESLAGLNLFLTQPVNGVPRVDFHWKFVTHPIASDNFYYPGYAGIGTQLQ
DSENELPPETTQQPNYESYSHRLSHIGLISASHVKALVYSWTHRSADRTNTIHSDSIT
QIPLVKAHTLQSGTTVVKGPGFTGGDILRRTSGGPFASFNVNLDWNLSQRYRARIRYA
STTNLRMYVTIAGERIFAGQFNKTMNTGDPLTFQSFSYATIDTAFTFPTKASSLTVGA
DTFSSGNEVYVDRFELIPVTATLEAVTDLERAQKAVHELFTSTNPGGLKTDVKDYHID
QVSNLVESLSDEFYLDEKRELFEIVKYAKQLHIEPNM"

gene 1..439
/gene="aroA"
/locus_tag="DR_1096"
Protein 1..439
/gene="aroA"
/locus_tag="DR_1096"
/product="**3-phosphoshikimate 1-carboxyvinyltransferase EPSPS**"
/EC_number="2.5.1.19«

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





pVS1 replicon for high stability in *Agrobacterium*

XhoI

Did you know this product can be customized?
find out more about customization options.

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





Close the microtubes and incubate at 37°C for 2h.


Reagents	HindIII/KpnI
10X buffer	4 μ l
DNA (1 μ g/ μ l)	4.0 μ l
<i>Hind</i> III (10u/ μ l)	2.0 μ l
KpnI	2.0 μ l
Water	28.0 μ l
Fin vol	40 μ l

<https://www.neb.com/tools-and-resources/interactive-tools/double-digest-finder?enzyme1={CB454480-8A3F-4BC8-843C-CB05AAFD0D57}&enzyme2={DA2537A1-428C-4E01-B72E-7E957F7A9654}>

Loading the Gel

Carefully place the pipette tip inside a well and gently/slowly expel the sample. The sample should sink into the well.

COMPONENT

20 µl REACTION

T4 DNA Ligase Buffer (10X)

2 µl

Vector DNA (4 kb)

50 ng (0.020 pmol)

Insert DNA (1 kb)

37.5 ng (0.060 pmol)

Nuclease-free water


to 20 µl

T4 DNA Ligase

1 µl

- Gently mix the reaction by pipetting up and down and microfuge briefly.
- For cohesive (sticky) ends, incubate at 16°C overnight or room temperature for 10 minutes.
- For blunt ends or single base overhangs, incubate at 16°C overnight or room temperature for 2 hours


Key Steps for Transformation


Plasmid DNA enters the bacterial cell and the genes are expressed.


- Bacterial cell suspension is placed in CaCl_2 solution
- Cells must be in *log phase* of growth.
- Cells are kept on ice until heat shock treatment
- Heat shock at 42 °C for one minute
- Recover period in LB broth
- Cells are spread on appropriate selection plates

CHEMICAL TRANSFORMATION WITH CALCIUM CHLORIDE

TRANSFORMATION BY ELECTROPORATION

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
Southern blot analysis of the selected transgenic events

Genomic DNAs prepared from GAI-10, 14, 56, and 71 plants were digested with *Bam*HI and then hybridized with a digoxigenin (DIG)-labelled probe specific to gene *G10*.

M: the DNA ladder


Fig. 3 Western blot analysis of transgene expression
Six plants were selected. Each sample was blotted with a primary antibody against Cry1Ac (a), Cry1Ig (b), and G10 (c), respectively. Prokaryotic expressed Cry1Ac, Cry1Ig, and G10 were used as a positive control (CK⁺), respectively. The sample prepared from non-transgenic rice was used as a negative control (CK⁻)




Fig. 4 Mortalities of cotton bollworm (CB), striped stem borer (SSB), and rice leaf roller (RLR) feeding on GAI-14

Non-transgenic rice at the same growing stage was used as the control (CK). Data are expressed as mean $\pm$ standard deviation ($n=3$). Mortalities of all insects feeding on GAI-14 plants were 100%. ** on the bars indicated extremely significant difference between CK and GAI-14 ($P<0.01$, Kruskal Wallis test)

The bioassay was conducted with cotton bollworm (a), striped stem borer (b), and rice leaf roller (c), respectively. Non-transgenic rice at the same growing stage was used as the control (CK)

The GAI-14 plants were sprayed with Roundup diluted at 1:140 (v/v). Non-transgenic rice at the same growing stage was used as a control (CK)

Table 2 Evaluation of basic agronomic performance between GAI-14 and Xiushui 134

Plant	Plant height (cm)	Panicles per plant	Grains per panicle	1000-grain weight (g)
GAI-14	69.05±2.24	11.43±2.49	112.53±6.16	24.55±0.67
Xiushui 134	68.31±1.97	12.10±2.53	110.15±8.32	24.22±0.49
<i>P</i> -value	0.19	0.34	0.40	0.50

The value was shown as mean±standard deviation ($n=30$). Data were analyzed by Student's *t*-test in SPSS

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ENZYMES USED IN MOLECULAR BIOLOGY


<u>Alkaline phosphatase</u>	Removes phosphate groups from 5' ends of DNA (prevents unwanted re-ligation of cut DNA)
<u>DNA ligase</u>	Joins compatible ends of DNA fragments (blunt/blunt or complementary cohesive ends). Uses ATP
<u>DNA polymerase I</u>	Synthesises DNA complementary to a DNA template in the 5'-to-3'direction. Starts from an oligonucleotide primer with a 3' OH end
Exonuclease III	Digests nucleotides progressively from a DNA strand in the 3' -to-5' direction
<u>Polynucleotide kinase</u>	Adds a phosphate group to the 5' end of double- or single-stranded DNA or RNA. Uses ATP
RNase A	Nuclease which digests RNA, not DNA
<u>Taq DNA polymerase</u>	Heat-stable DNA polymerase isolated from a thermostable microbe (<i>Thermus aquaticus</i>)

RESTRICTION ENZYMES

- **The restriction enzymes most used in molecular biology labs cut within their recognition sites and generate one of three different types of ends.**


5' OVERHANGS

- **5' overhangs:** The enzyme cuts asymmetrically within the recognition site such that a short single-stranded segment extends from the 5' ends. *Bam* HI cuts in this manner.


3' OVERHANGS

- **3' overhangs:** Again, we see asymmetrical cutting within the recognition site, but the result is a single-stranded overhang from the two 3' ends. *KpnI* cuts in this manner.


BLUNT ENDS

- **Blunts:** Enzymes that cut at precisely opposite sites in the two strands of DNA generate blunt ends without overhangs. *Sma*I is an example of an enzyme that generates blunt ends.


Converting a 5' overhang to blunt end

- Both Klenow and T4 DNA polymerase can be used to fill in 5' protruding ends with dNTPs
- Used in joining DNA fragments with incompatible ends
- Once the ends have been blunted, ligation can proceed

Converting a 3' overhang to a blunt end

- T4 DNA polymerase has a 3'-5' exonuclease activity
- In the presence of excess dNTPs will convert a 3' protruding end to a blunt end
- Ligation can now proceed


DIRECTIONAL CLONING

- Often one desires to insert foreign DNA in a particular orientation
- This can be done by making two cleavages with two different **restriction** enzymes
- Construct foreign DNA with same two **restriction** enzymes
- Foreign DNA can only be inserted in one direction




- Good efficiency of ligation of foreign DNA into a vector can be achieved if both the vector and the insert DNA are cut with 2 different restriction enzymes which leave single stranded ends (cohesive ends).
- The DNA is ligated in only one direction, and there is only a low background of non-recombinant plasmids.
- If only one restriction enzyme is used to cut the vector and insert, then efficiency of ligation is lower, DNA can be inserted in two directions and tandem copies of inserts may be found.
- To avoid high background of non-recombinants, **alkaline phosphatase** is used to remove 5' phosphate groups from the cut vector to prevent self-ligation.

Alkaline phosphatase

- Alkaline phosphatase removes 5' phosphate groups from DNA and RNA. It will also remove phosphates from nucleotides and proteins. These enzymes are most active at alkaline pH


ENZYMES USED IN MOLECULAR BIOLOGY


Alkaline phosphatase

- There are two primary uses for alkaline phosphatase in DNA manipulations:
- **Removing 5' phosphates from plasmid and bacteriophage vectors** that have been cut with a restriction enzyme. In subsequent ligation reactions, this treatment prevents self-ligation of the vector and thereby greatly facilitates ligation of other DNA fragments into the vector (e.g. subcloning).
- **Removing 5' phosphates from fragments of DNA prior to labeling with radioactive phosphate.** Polynucleotide kinase is much more effective in phosphorylating DNA if the 5' phosphate has previously been removed

DEPHOSPHORYLATED VECTOR


R.E.S WITH COMPATIBLE ENDS


Generating a new R.E. site at a blunt end

- **Use linkers to generate a new R.E.**
- **Linkers are used to place sticky ends on to a blunt-ended molecule**
- **Short blunt ended synthetic ds DNA containing a R.E. site**
- **Experimental design:**
 - (i) Blunt ended DNA + linker (T4 ligase)**
 - (ii) Digest with appropriate R.E.**
 - (iii) Ligate to vector**



Introducing a R.E. site by PCR



- R.E. site is designed into the 5' end of the PCR primer
- PCR fragment is digested with appropriate R.E., purified and ligated into plasmid vector

USING DIFFERENT R.E.s

Anneal and extend oligos with 5' restriction sites
(EcoRI on one end, BamHI on the other)

Correct
Orientation

Wrong
Orientation